



SOL-CARE™ Safety Needles and Safety Needle Combinations



Features & Benefits

- Easy to use, easy to handle while wearing gloves
- One-handed technique for activation of the integrated safety shield via 3 methods: finger; thumb and hand; or hard surface
- The color-coded safety shield and needle hub provide easy identification of the needle size
- The bevel up needle allows for positive identification of the bevel position and facilitates low angled injections without interference of the safety cover
- Both hands remain behind the needle during activation of the safety feature on hard surface
- An audible click and visible change indicate activation of the safety cover
- The needle remains securely locked with a double-lock mechanism once activated to provide added security
- Compatible with both standard luer lock and luer slip tip syringes



SOL-CARE™ Safety Needle



Ref	Gauge	Length	Qty/Box	Qty/Case
SN3012	30	1/2"	100	800
SN27125	27	1 1/4"	100	800
SN2758	27	5/8"	100	800
SN2712	27	1/2"	100	800
SN2558	25	5/8"	100	800
SN2510	25	1"	100	800
SN2515	25	1 1/2"	100	800
SN2315	23	1 1/2"	100	800
SN2310	23	1"	100	800
SN2210	22	1"	100	800
SN2215	22	1 1/2"	100	800
SN2110	21	1"	100	800
SN2115	21	1 1/2"	100	800
SN2010	20	1"	100	800
SN2015	20	1 1/2"	100	800
SN1915	19	1 1/2"	100	800
SN1810	18	1"	100	800
SN1815	18	1 1/2"	100	800



SOL-CARE™ Safety Needle Combinations



Ref	Size	Gauge	Length	Qty/Box	Qty/Case
32558SN	3mL	25	5/8"	100	1200
32510SN	3mL	25	1"	100	1200
32310SN	3mL	23	1"	100	1200
32210SN	3mL	22	1"	100	1200
32215SN	3mL	22	1 1/2"	100	1200
32110SN	3mL	21	1"	100	1200
32115SN	3mL	21	1 1/2"	100	1200
32010SN	3mL	20	1"	100	1200
32015SN	3mL	20	1 1/2"	100	1200